

## SpaceSAVER Bouncy Climber™



### Introduction

Designed to withstand heavy, clinical use, the SpaceSAVER Bouncy Climber™ is a child AND adult favorite. It's durable construction makes it safe to use as a bouncer for a child or adult up to 175 lbs, and it's shape lets it fit easily into the corner of a room.

The yellow bounce-surround keeps bouncers centered and safe, giving piece of mind—especially when working with children with poor impulse control. Have kids that need calming? The rhythmic, calming, up and down movement with heavy impact to receptors in the body yields an organized child that is ready to learn!

### Product Use


**WARNING:** The SpaceSAVER Bouncy Climber is intended to be used under the supervision of an adult at all times. Depending on the activity undertaken, proper floor matting and spotting techniques should be used. As with any piece of equipment, failure to attend to a child during use may result in injury.

**WARNING:** This is an inflatable product. Take care to keep sharp objects out of contact with it, which will void your warranty. If used with shoes, inspect shoes for foreign objects or sharp elements prior to use.


**To inflate:** The SpaceSAVER Bouncy Climber™ has two valves, one on the bottom (blue material), and one on the top (yellow material). These valves are spring loaded so that they will close automatically to keep the air in as you are filling it.



To open, rotate the valve counter-clockwise and pull



Set the outer plug to the side. Note the spring loaded plunger in the center.



Press your inflator's hose onto the valve and inflate. (See recommended products)



To deflate, press plunger in. To hold valve open, press plunger and rotate clockwise.

**Note:** If you use the hold-open feature of the valve to deflate, make sure to close the valve before re-inflating. Otherwise, air will escape when re-filling, making it difficult to fully inflate. To close the valve, fully depress the center plunger and rotate it counter-clockwise. The plunger will pop up and the valve will close.

(continued...)

### To use as a Bouncer:

- Arrange the SpaceSAVER Bouncy Climber™ where appropriate for your space, blue side down (see image, previous page). You will notice that one side of the bouncer does not have a lower yellow “rib”. This is the open side and is the easiest way to get onto the bouncer.


### To use as a “nest”:

- Turn the SpaceSAVER Bouncy Climber™ yellow side down. The bouncing surface is now the “floor” of the nest (see image, right).  
**Warning: The center of gravity is higher when arranged in this manner. Supervise carefully when used in this configuration. We recommend using the provided anchor straps to secure the “nest” to the wall (see image “A”, below).**



### To anchor a “nest” to the wall:

- Place the SpaceSAVER Bouncy Climber™ in a location against the wall where it will be used. Note the sets of “D” rings on the blue material. Mark the wall next to both sets of “D” rings. Using appropriate anchors for your type of wall (brick / wood / drywall, not included), attach the provided webbing straps to the wall through the grommet (use large washers to make sure the webbing does not pull through the fastener, image B). Next, thread the webbing straps through the “D” rings as shown (image C-E). **Note: This provides a good deal of stability when arranged as a “nest”, but does not eliminate the need for proper supervision and/or perimeter floor mats.**



### To attach the Wonder Wedge to the Bouncer:

- The Wonder Wedge™ (see recommended Products) has webbing straps that were designed to connect with the SpaceSAVER Bouncy Climber™. Simply lean the Wonder Wedge™—straps facing down—against any side of the Bouncer with “D” rings, and connect in the same manner as illustrated above.

## Safety

- Even though the SpaceSAVER Bouncy Climber™ is one of the safest trampoline-like products around, it must ALWAYS be used with supervision. Depending on the needs and capabilities of the child and your assessment as a therapist, teacher, or parent, proper protection and safety devices—such as a helmet—may be required.
- When used inverted as a “nest”, use adequate perimeter matting and proper spotting techniques to ensure the safety of the child.
- In general, the safety of any therapeutic situation is influenced by matching appropriate physical challenges to the needs and abilities of the child, as well as the layout, condition, and maintenance of your room and equipment.
- These are general safety guidelines. In addition, please refer to the safety policies of your facility, organization, and/or state, as these guidelines are not intended to supersede or replace them.
- Regular inspection is a key part of any safety program. Please utilize the included inspection checklist to ensure best safety practices.

## Cleaning and Storage

*Please reference your facility's standards and regulations on cleaning and disease prevention before using these guidelines.*

- The SpaceSAVER Bouncy Climber™ can be cleaned with a mild soap and water solution and/or disinfectant.
- Dry before using or storage.
- Store the SpaceSAVER Bouncy Climber™ in a dry, climate controlled space, out of direct sunlight. It is preferable to keep this product out of extreme heat for storage.

## Recommended Products

Wonder Wedge: Item # 1359001

Heavy Duty 2.5 HP Inflator: Item # 006785

Super Heavy Duty Inflator/Deflator: Item # 003743

Fat mat: Item # 013434

Space Walk: Item # 031551

Striped Super Expand-O-Mats: many sizes available—see website



©2010 School Specialty, Inc. • Greenville, WI • USA • All rights reserved  
Customer Service **1-888-388-3224**

Or reach us on the web at: <http://www.schoolspecialty.com>

**1-YEAR GUARANTEE OF SATISFACTION**

*See more fun activity guides for other great Abilitations products at [www.schoolspecialty.com](http://www.schoolspecialty.com)*

